

Product datasheet for **MR226858**

Lepr (NM_001122899) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lepr (NM_001122899) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lepr
Synonyms:	db; diabetes; Leprb; LEPROT; Modb1; OB-RGRP; obese-like; obl; Obr
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR226858 representing NM_001122899
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGTGTCAGAAATTCTATGTGGTTTTGTACTGGAATTTCTTTATGTGATAGTGCACCTTAACC
 TGGCATATCCAATCTCTCCCTGGAAATTTAAGTTGTTTTGTGGACCACCGAACACAACCGATGACTCCTT
 TCTCTCACCTGTGGAGCCCCAAACAATGCCTCGGCTTTGAAGGGGGCTTCTGAAGCAATTGTTGAAGCT
 AAATTTAATCAAGTGGTATCTACGTTCTGAGTTATCCAAAACAGTCTTCCACTGTTGCTTTGGGAATG
 AGCAAGGTCAAACCTGCTCTGCACTCACAGACAACACTGAAGGGAAGACACTGGCTTCAGTAGTGAAGGC
 TTCAGTTTTTCGCCAGCTAGGTGTAACACTGGGACATAGAGTGTGGATGAAAGGGGACTTGACATTATTC
 ATCTGTCATATGGAGCCATTACCTAAGAACCCTTCAAGAATTATGACTCTAAGGTCCATCTTTTATATG
 ATCTGCCTGAAGTCATAGATGATTCGCCTCTGCCCCACTGAAAGACAGCTTTCAGACTGTCCAATGCAA
 CTGCAGTCTTCGGGGATGTGAATGTCATGTGCCGGTACCCAGAGCCAACTCAACTACGCTCTTCTGATG
 TATTTGGAATACATCTGCCGGTGTGAGTTTTAGTCACCTCTGATGTCAGCTGCAGCCATGCTTGTG
 TGAACCCCGATCCACCCTTAGGTTTGCATATGGAAGTCACAGATGATGGTAAATTTAAAGATTTCTTGGGA
 CAGCCAAACAATGGCACCATTTCCGCTTCAATATCAGGTGAAATTTAGAGAATTCTACAATTGTAAGA
 GAGGCTGCTGAAATGTCTCAGCTACATCTCTGCTGGTAGACAGTGTGCTTCTGGATCTTCATATGAGG
 TCCAGGTGAGGAGCAAGAGACTGGATGGTTCAGGAGTCTGGAGTACTGGAGTTCACCTCAAGTCTTTAC
 CACACAAGATGTTGTGTATTTCCACCCAAAATCTGACTAGTGTGGATCGAATGCTTCTTTTTCATTGC
 ATCTACAAAAACGAAAACCAGATTATCTCTCAAACAGATAGTTGGTGGAGGAATCTAGCTGAGAAAA
 TCCCTGAGATACAGTACAGCATTGTGAGTGACCGAGTTAGCAAAGTTACCTTCTCAACCTGAAAGCCAC
 CAGACCTCGAGGGAAGTTTACCTATGACGCAAGTACTGCTGCAATGAGCAGGCGTCCATCACCCTAT
 GCTGAATTATACGTGATCGATGTCAATATCAATATATCATGTGAAACTGACGGGACTTAACTAAAATGA
 CTTGCAGATGGTCACCCAGCACAATCCAATCACTAGTGGGAAGCACTGTGCAGCTGAGGTATCACAGGCG
 CAGCCTGTATTGCTGTAGTCCATCTATTATCCTACGTCTGAGCCAAAAACTGCGTCTTACAGAGA
 GACGGCTTTTATGAATGTGTTTTCCAGCCAATCTTTCTATTATCTGGCTATACAATGTGGATCAGGATCA
 ACCATTCTTTAGGTTCACTTGACTCGCCACCAACGTGTCTTCTGACTCCGTAGTAAAACCACTACC
 TCCATCTAACGTAAAAGCAGAGATTACTGTAACACTGGATTATTGAAAGTATCTTGGGAAAAGCCAGTC
 TTTCCGGAGAATAACCTTCAATCCAGATTCGATATGGCTTAAGTGGAAAAGAAATACAATGGAAGACAC
 ATGAGGTATTCGATGCAAAGTCAAAGTCTGCCAGCCTGCTGGTGTGAGACCTCTGTGCACTATGTGGT
 CCAGGTTTCGCTGCCGGCGGTTGGATGGACTAGGATATTGGAGTAATTGGAGCAGTCCAGCCTATACGCTT
 GTCATGGATGTAAGGTTCTATGAGAGGGCTGAATTTTGGAGAAAAATGGATGGGGACGTTACTAAAA
 AGGAGAGAAAATGTCACCTTGCTTTGGAAGCCCCTGACGAAAAATGACTCACTGTGTAGTGTGAGGAGGTA
 CGTGGTGAAGCATCGTACTGCCACAATGGGACGTGGTCAGAAGATGTGGGAAATCGGACCAATCTCACT
 TTCCTGTGGACAGAACCAGCGCACACTGTTACAGTTCTGGCTGTCAATCCCTCGGCGTTCCTTGTGA
 ATTTTAACTTACCTTCTCATGGCCATGAGTAAAGTGTGCTGTGGAGTCACTCAGTGCTTATCCCT
 GAGCAGCAGCTGTGTCATCCTTTCTGGACACTGTACCTGATGATTATAGTCTGTTATATCTGGTTATT
 GAATGGAAGATCCTTAATGAAGATGATGGAATGAAGTGGTTAGAATCCCTCGAATGTTAAAAAGTTTT
 ATATCCACGATAATTTTATCCCATCGAATAATCAGTTTAGTCTTTACCCAGTATTTATGGAAGGAGT
 TGGAAAACCAAAGATAATTAATGGTTTACCAAAGATGCTATCGACAAGCAGCAGAATGACGCAGGGCTG
 TATGTCATTGTACCCATAATTTTCTTGTGCTCTACTGCTCGGAACACTGTTAATTTACACCAGA
 GAATGAAAAAGTTGTTTTGGGACGATGTTCAAACCCCAAGAATTGTTCTGGGCACAAGGACTGAATTT
 CCAAAGAGAACGGACTCTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR226858 representing NM_001122899
Red=Cloning site Green=Tags(s)

```

MMCQKFYVLLHWEFLYVIAALNLAYPISPWKFKLCGPPNTTDDSFLSPAGAPNNASALKGASEAIVEA
KFNSSGIYVPELSTVVFHCCFGNEQQNCSTLDNTEGKTLASVVKASVFRQLGVNWDIECWMKGDTLF
ICHMEPLPKNPFKNYDSKVHLLYDLPEVIDDSPLPLKDSFQTVQCNCSLRGCECHVPVPRAKLNYALLM
YLEITSAGVSVFQSPMLSLQPMPLVVKPDPPLGLHMEVTDGDLKISWDSQTMAPPFLQYQVKYLENSTIVR
EAAEIVSATSLLVDSVLPGSSYEVQVRSKRLDGSVWSDWSSPQVFTTQDVVYFPFKILT SVGSNASFHC
IYKNENQIISSEKQIVVWRNLAEKIPEIQYSIVSDRVSKVTFSNLKATRPRGKFTYDAVYCCNEQACHRY
AELYVIDVNIINISCETDGYLTKMTCRWSPSTIQSLVGSTVQLRYHRRSLYCPDPSIHTPSEPKNVQLQR
DGFYECVFPQIFLLSGYTMWIRINHSLGSLDSPPTCVLPDSVVKPLPPSNVKAETVNTGLLKVSWEKPV
FPENNLQFQIRYGLSGKEIQWKTHEVFDKSKSASLLVSDLCVAVVVQVRCRRLDGLGYWSNWSSPAYTL
VMDVKVPMRGPWFWRKMDGDVTKKERNVTLWKPLTKNDSLCSVRRYVVKHRTAHNGTSEDVGNRNTLT
FLWTEPAHTVTVLAVNSLGASLVNFNLTFSWPMSKVSAYESLSAYPLSSSCVILSWTLSPDDYSLLYLVI
EWKILNEDDGMKWLRIIPSNVKKFYIHDNFIPIEKYQFSLYPVFMGVGPKIINGFTKDAIDKQONDAGL
YVIVPIIISSEVLLLGLLLISHQRMKCLFWDDVNPKNCSWAQGLNFQKRTDTL
    
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9003_d02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF


ACCN: NM_001122899

ORF Size:	2682 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001122899.2
RefSeq Size:	5571 bp
RefSeq ORF:	2685 bp
Locus ID:	16847
UniProt ID:	P48356
Cytogenetics:	4 46.96 cM
MW:	101.5 kDa

Gene Summary:

Receptor for hormone LEP/leptin (Probable) (PubMed:11861497). On ligand binding, mediates LEP central and peripheral effects through the activation of different signaling pathways such as JAK2/STAT3 and MAPK cascade/FOS (PubMed:10799542, PubMed:25383904, PubMed:11923481, PubMed:11861497). In the hypothalamus, LEP acts as an appetite-regulating factor that induces a decrease in food intake and an increase in energy consumption by inducing anorexigenic factors and suppressing orexigenic neuropeptides, also regulates bone mass and secretion of hypothalamo-pituitary-adrenal hormones (PubMed:10660043, PubMed:12594516). In the periphery, increases basal metabolism, influences reproductive function, regulates pancreatic beta-cell function and insulin secretion, is pro-angiogenic and affects innate and adaptive immunity (PubMed:25383904, PubMed:11923481). Control of energy homeostasis and melanocortin production (stimulation of POMC and full repression of AgRP transcription) is mediated by STAT3 signaling, whereas distinct signals regulate NPY and the control of fertility, growth and glucose homeostasis (PubMed:12594516). Involved in the regulation of counter-regulatory response to hypoglycemia by inhibiting neurons of the parabrachial nucleus (PubMed:25383904). Has a specific effect on T lymphocyte responses, differentially regulating the proliferation of naive and memory T-cells. Leptin increases Th1 and suppresses Th2 cytokine production (PubMed:9732873).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR226858