

Product datasheet for **MR226863**

Olfr16 (NM_008763) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr16 (NM_008763) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Olfr16
Synonyms:	MOR23; MOR26713
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR226863 representing NM_008763 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCAGAGAAATAACTTCACTGAAGTGATAGAGTTCGTCTTCCTGGGATTCTCCAGCTTTGGAAAGCATC
AGATAACCCTCTTTGTGGTTTTCTAACCATCTACATTTAACTCTGGCTGGCAACATCATTATAGTGAC
AATCACACACATAGACCACCACCTTCACACTCCCATGTACTTCTTTCTGAGCATGTTGGCAAGCTCAGAG
ACTGTGTACACACTGGTCATTGTCCACGGATGCTTTCCAGCCTGATTTTTTACAACCTCCCATATCCT
TGGCAGGCTGCGCAACCCAAATGTTCTTTTTGTACCTTGGCCACCAACAAGTCTTTCTGCTCACAGC
AATGGGTTATGATCGTTATGTGGCTATTTGTAATCCTCTGAGATATACAATCATCATGAGCAAGGGAATG
TGTGCCTTGTTGGTTTGTGGGTCTTTAGGCACTGGCCTGGTTATGGCAGTTCCTCATGTGCCAGCCATGT
TCCATTTGCCCTTTTGTGGCACGGTGGTGGAGCACTTTTTCTGTGACATATACCCAGTAATGAAGCTTTC
TTGTGTTGATACCACTGTCAATGAGATAATCAATTATGGTGTAAGTTCATTTGTAATTCCTGTGCCATA
GGGCTGATTTTATCTCCTATGTGCTCATTGTCTCTTCCATCCTTAAATTTGTGTCCACTGAAGGCCAGA
AGAAAGCCTTTGCCACCTGTGCCTCTCATCTCACTGTGGTCATTGTCCACTATGGCTGTGCCTCCATTGC
CTACCTCAAACCCAAATCAGAAAGTTCAAGTAAAAAGACCTTCTCTCTGTGACCTACACCATCATC
ACTCCCTTGCTGAACCCTGTTGTCTACAGCCTCAGGAACAAAGAAGTCAAAGATGCTCTATGCAGAGCTG
TGGGCAGAAACACTTCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR226863 representing NM_008763
Red=Cloning site Green=Tags(s)

MQRNNFTEVIEFVFLGFSFSGKHQITLFFVFLTIYILTLAGNIIIVTITHIDHHLHTPMYFFLSMLASSE
TVYTLVIVPRMLSSLIFYNLPISLAGCATQMFVFVTLATNNCFLTAMGYDRYVAICNPLRYTIIMSKGM
CALLVCGSLGTGLVMAVLHVPAMFHLPCGTVVEHFFCDIYPVMKLSVDTTVNEIINYGVSSFVILVPI
GLIFISYVLIVSSLLIKVSTEGQKAFATCASHLTVVIVHYGCASIAYLKPKSESSVEKDLLLSVTYTI
I
TPLLNPVYVSLRNKEVKDALCRAVGRNTS



TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mm9023_b04.zip

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_008763

ORF Size: 927 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_008763.2](#), [NP_032789.1](#)

RefSeq Size: 1050 bp

RefSeq ORF: 930 bp


Locus ID: 18313

Cytogenetics: 1 H3

MW: 34.7 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for MR226863