

Product datasheet for **MR226870**

Cry2 (NM_001113333) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cry2 (NM_001113333) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cry2
Synonyms:	AV006279; D130054K12Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR226870 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCGCGGCTGCTGTGGTGGCAGCGACGGTCCCCGCGCAATCGATGGGCGCGGACGGCGCGTCTCTCGG
 TGCACTGGTTCGCAAAGGACTACGGCTCCACGACAACCCGCGCTGCTAGCTGCCGTGCGCGGGCGCG
 CTGTGTGCGCTGCGTCTACATCCTCGACCCGTGGTTCGCGGCTCCTCGTCTGTGGGCATCAACCGATGG
 AGGTTCTACTGCAATCTCTGGAAGATCTGGACACAAGCTTAAGAAAGCTGAATTCGCGTCTGTTTGTAG
 TCCGGGGACAGCCAGCTGATGTGTTCCCAAGGCTGTTCAAGGAATGGGGGTGACCCGTTTGACCTTTGA
 ATATGACTCTGAACCTTTGGGAAAGAACGGGATGCAGCCATTATGAAGATGGCCAAGGAGGCTGGCGTG
 GAGGTGGTACTGAGAACTCTCACACCCTATGACCTAGACAGAATCATCGAAGTGAATGGGCAGAAAC
 CACCCCTTACCTACAAGCGCTTTCAGGCCCTCATCAGCCGATGGAGCTGCCAAGAAGCCCGGGTGGC
 TGTGAGCAGCCAGAGATGGAGAGCTGCAGAGCTGAGATCCAGGAGAACCATGACGACACCTATGGCGTG
 CCTTCCCTGGAGGAGCTGGGATTCACCCAGGAAGGACTTGGCCAGCTGTTTGGCAAGGAGGAGAGACAG
 AAGCTCTGGCCCGCTGGACAAGCACTTGAACCGAAGGCTGGGTTGCCAACTATGAGAGACCTCGGAT
 GAATGCCAATTCCTACTGGCCAGCCCCACAGGCTCAGCCCTACCTGCGCTTTGGATGCCTCTCTCTGC
 CGCCTCTTCTACTACCGCTGTGGGACTTGTACAAGAAGTGAAGAGGAACAGCACACCCCCCTCTCTCT
 TATTTGGACAACCTCTGTGGCGAGAATTCTTCTACACAGCGGCCACCAACAACCCAGGTTTGACCGAAT
 GGAGGGAAACCCATCTGCATCCAGATCCCCTGGGACCGCAACCCGAAGCCCTGGCCAAGTGGGCCGAG
 GGCAAGACAGGCTTCCCTTGGATTGACGCCATCATGACCAACTGAGGCAGGAGGGCTGGATCCACCACC
 TGGCCCGGACGCTGTGGCTGCTTCTCACCCGCGGGACCTCTGGGTGAGTGGGAGCGGGGCTCCG
 GGATTTGACGAGCTGCTCTGGATCCGATTTTCAAGTGTGAATGCAGGCAGCTGGATGTGGCTGTCCTGC
 AGTGCTTTCTTCCAACAATTCTTCCACTGCTACTGCCCTGTGGGCTTCGGCCGACGTACAGACCCCAAGT
 GGGACTACATCCGGCGGTACCTGCCAACTGAAAGGCTTCCCCTCTCGATACATCTATGAGCCCTGGAA
 TGCCCCGAGTCAGTTCAGAAGGCTGCCAAGTGCATCATTGGCGTGGACTACCCACGGCCCATCGTCAAT
 CATGCAGAGACTAGTCGGCTCAACATTGAACGAATGAAGCAGATCTACCAACAGCTGTCGAGATACCGGG
 GACTCTGTCTATTGGCATCTGTCCCTTCTGTGTGAAGACCTCAGTACCCTGTGGCAGAGCCTGGTTC
 AAGCCAAGCTGGGAGCATCAGCAACACAGGCCCCAGAGCACTATCCAGTGGCCAGCTTCCCCAAACGC
 AAGCTGGAAGCAGCCGAGGAACCTCTGGTGAAGAAGTACCAAGCGGGCTAGAGTGACGGAGATGCCTA
 CCAAGAGCCAGCAAGCAAGGACTCC

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR226870 protein sequence
 Red=Cloning site Green=Tags(s)

MAAAVVAATVPAQSMGADGASSVHFRKGLRLHDPALLAAVARGARCVRCVYILDPAFASSSVGINRW
 RFLQLSLEDLDTSLRKLNSRLFVVRGQPADVFPRLFKEWGVTRLTFEYDSEPFGKERDAAIMKMAKEAGV
 EVVTENSHLYDLDRIIELNGQKPLTYKRFQALISRMELPKKPAVAVSSQOMESCRAEIQENHDDTYGV
 PSLEELGFPTGLGPAVWQGGETEALARDKHLERKAWVANYERPRMNANSLASPTGLSPYLRFGCLSC
 RLFYYRLWDLYKVKRNSTPPLSLFGQLLWREFFYTAATNNPRFDRMEGNPICIQIPWDRNPEALAKWAE
 GKTGFPPWIDAIMTQLRQEGWIHHLARHAVACFLTRGDLWVSWESGVRVFDLLEDADFVNAGSMMWLS
 SAFFQQFFHCYCPVGFGRRTDPSGDYIRRYLPLKKGFPSTRYIYEPWNAPESVQKAAKCIIGVDYPRPIVN
 HAETSRLNIERMKQIYQQLSRYRGLCLLASVPSCVEDLSHPVAEPGSSQAGSISNTGPRALLSSGPASP
 KLEAAEPPGEELTKRARVTEMPTQEPASKDS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001113333

ORF Size: 1779 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001113333.1](#), [NP_001106804.1](#)

RefSeq Size: 2193 bp

RefSeq ORF: 1779 bp


Locus ID: 12953

Cytogenetics: 2 E1

MW: 66.9 kDa

Gene Summary: This gene encodes a flavin adenine dinucleotide-binding protein that is a key component of the circadian core oscillator complex, which regulates the circadian clock. This gene is upregulated by Clock/Arntl heterodimers but then represses this upregulation in a feedback loop using Per/Cry heterodimers to interact with Clock/Arntl. Polymorphisms in this gene have been associated with altered sleep patterns. The encoded protein is widely conserved across plants and animals. [provided by RefSeq, Feb 2014]

Product images:



Circular map for MR226870