

Product datasheet for **MR226873L3V**

Tap1 (NM_013683) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tap1 (NM_013683) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tap1
Synonyms:	ABC17; Abcb; Abcb2; APT1; Ham; Ham-; Ham-1; Ham1; MTP; MTP1; PSF; PSF1; RI; RING4; T; TAP; Tap-1; Y3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013683
ORF Size:	2172 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226873).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013683.2 , NP_038711.2
RefSeq Size:	2950 bp
RefSeq ORF:	2175 bp
Locus ID:	21354
UniProt ID:	P21958
Cytogenetics:	17 17.98 cM



[View online »](#)

Gene Summary:

The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. This protein forms a heterodimer with Tap2 that transports short peptides from the cytosol into the endoplasmic reticulum lumen. Mutations in the human gene may be associated with ankylosing spondylitis, insulin-dependent diabetes mellitus, and celiac disease. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2009]