

## Product datasheet for **MR226881L4V**

### Socs2 (NM\_001168657) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Socs2 (NM_001168657) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Socs2
Synonyms:	8030460M17; AI527257; AW108012; CIS2; Cish2; D130043N08Rik; hg; JAB; SOCS-2; SSI-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001168657
ORF Size:	594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226881).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001168657.1</a> , <a href="#">NP_001162128.1</a>
RefSeq Size:	2074 bp
RefSeq ORF:	597 bp
Locus ID:	216233
UniProt ID:	<a href="#">O35717</a>
Cytogenetics:	10 49.35 cM



[View online »](#)

**Gene Summary:**

SOCS family proteins form part of a classical negative feedback system that regulates cytokine signal transduction. SOCS2 appears to be a negative regulator in the growth hormone/IGF1 signaling pathway. Probable substrate recognition component of a SCF-like ECS (Elongin BC-CUL2/5-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins (By similarity). [UniProtKB/Swiss-Prot Function]