

## Product datasheet for **MR226896L4V**

### **Ctla4 (NM\_009843) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ctla4 (NM_009843) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ctla4
Synonyms:	Cd152; Ctla; Ctla-4; Ly-5; Ly-56
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_009843
ORF Size:	669 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226896).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_009843.4</a> , <a href="#">NP_033973.2</a>
RefSeq Size:	1933 bp
RefSeq ORF:	672 bp
Locus ID:	12477
UniProt ID:	<a href="#">P09793</a>
Cytogenetics:	1 30.58 cM



[View online »](#)

**Gene Summary:**

This gene is a member of the immunoglobulin superfamily, and encodes a protein that functions as a negative regulator of T-cell responses. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Aug 2013]