

Product datasheet for **MR226905L3V**

Csf3r (NM_007782) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Csf3r (NM_007782) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Csf3r
Synonyms:	Cd114; Csfgr; G-CSF-R; G-CSFR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_007782
ORF Size:	2514 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226905).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007782.2
RefSeq Size:	3907 bp
RefSeq ORF:	2514 bp
Locus ID:	12986
UniProt ID:	P40223
Cytogenetics:	4 59.77 cM
Gene Summary:	Receptor for granulocyte colony-stimulating factor (CSF3). In addition it may function in some adhesion or recognition events at the cell surface.[UniProtKB/Swiss-Prot Function]



[View online »](#)