

Product datasheet for **MR226942A1V**

Mouse Defa2 (NM_001195634) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Defa2 (NM_001195634) AAV Particle
Tag: Myc-DDK
Symbol: Defa2
Synonyms: CR2; Defcr-2; Defcr2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR226942 representing NM_001195634
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAAGACACTTGTCTCTCTCTGCCCTCGTCTGCTGGCCTCCAGGTCCAGGCTGATCCTATCCAAA
ACACAGATGAGGAGACTAATACTGAGGAGCAGGCAGGGGAGGAGGACCAGGCTGTGTCTGTCTCTTTGG
AGACCCAGAAGGATCTGCTCTTCATGAAAAATCGTCGAGAGATCTGATATGCTATTGTAGAAAAGGAGGC
TGCAATAGAGGAGAACAAGTTTATGGACCTGCTCAGGACGACTTTTGTCTGCTGCCGCCGCCGCCACC
GCCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAAGTTTAA

Protein Sequence: >MR226942 representing NM_001195634
Red=Cloning site **Green**=Tags(s)

MKTLVLLSALVLLAFQVQADPIQNTDEETNTEEQAGEEDQAVSVSFGDPEGSALHEKSSRDILICYRKGK
CNRGEQVYGTCSGRLLFCCRHRH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse
Serotype: AAV-2
ACCN: NM_001195634
ORF Size: 285 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001195634.2, NP_001182563.2</u>
RefSeq Size:	424 bp
RefSeq ORF:	288 bp
Locus ID:	100294660
Cytogenetics:	8
MW:	10.6 kDa