

Product datasheet for **MR226998L3V**

Nodal (NM_013611) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nodal (NM_013611) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nodal
Synonyms:	Tg.413d
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013611
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR226998).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013611.4
RefSeq Size:	2094 bp
RefSeq ORF:	1065 bp
Locus ID:	18119
UniProt ID:	P43021
Cytogenetics:	10 32.21 cM



[View online »](#)

Gene Summary:

This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate the mature protein, which regulates early embryonic development. Homozygous knockout mice for this gene exhibit early embryonic lethality, while expression of a hypomorphic allele results in defects in anteroposterior and left-right patterning. [provided by RefSeq, Aug 2016]