

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for MR227022L4V

Gata4 (NM_008092) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gata4 (NM_008092) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gata4
Synonyms:	Gata-4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008092
ORF Size:	1326 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227022).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 008092.3, NP 032118.2</u>
RefSeq Size:	3393 bp
RefSeq ORF:	1326 bp
Locus ID:	14463
UniProt ID:	<u>Q08369</u>
Cytogenetics:	14 33.24 cM




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Gata4 (NM_008092) Mouse Tagged ORF Clone Lentiviral Particle – MR227022L4V

Gene Summary:Transcriptional activator that binds to the consensus sequence 5'-AGATAG-3' and plays a key
role in cardiac development (By similarity). In cooperation with TBX5, it binds to cardiac
super-enhancers and promotes cardiomyocyte gene expression, while it downregulates
endocardial and endothelial gene expression (By similarity). Involved in bone morphogenetic
protein (BMP)-mediated induction of cardiac-specific gene expression (By similarity). Binds to
BMP response element (BMPRE) DNA sequences within cardiac activating regions (By
similarity). Acts as a transcriptional activator of ANF in cooperation with NKX2-5
(PubMed:9584153). Promotes cardiac myocyte enlargement (By similarity). Required during
testicular development (By similarity). May play a role in sphingolipid signaling by regulating
the expression of sphingosine-1-phosphate degrading enzyme, spingosine-1-phosphate lyase
(By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



[MR227022L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR227022L4V particle to overexpress human Gata4-mGFP fusion protein.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US