

## Product datasheet for **MR227075L4V**

### Elavl4 (NM\_001038698) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Elavl4 (NM_001038698) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Elavl4
Synonyms:	Elav; Hud; PNEM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001038698
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227075).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001038698.1</a> , <a href="#">NP_001033787.1</a>
RefSeq Size:	4008 bp
RefSeq ORF:	1143 bp
Locus ID:	15572
UniProt ID:	<a href="#">Q61701</a>
Cytogenetics:	4 51.42 cM
Gene Summary:	May play a role in neuron-specific RNA processing. Protects CDKN1A mRNA from decay by binding to its 3' UTR. Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF and FOS mRNA (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)