

Product datasheet for MR227076

Rpl29 (NM_009082) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rpl29 (NM_009082) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rpl29
Synonyms: Rpl43
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR227076 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCAAGTCCAAGAACCACACCACACAACCAGTCCCCAAATGGCACAGAAATGGCATCAAGAAAC
CCCGGTCGCAAAGATACGAATCTCTTAAGGGGTTGACCCCAAGTTCCTGAGGAACATGCGCTTTGCCAA
GAAGCACAAAGAAGGCTGAAGAAGATGCAGGCCAACATGCAAAGGCAGTGAAGTGCAGCGCAGAG
GCCATCAAGGCCCTGGTGAAGCCTCAGGCCATCAAGCCAAAGATGCCAAAAGGCCCAAACCTCAAGCGG
TGGCTTTCATCGCTCACCCAAGCTTGGGAAGCGGATTCGAAGCTACATGGCCAAGGGTCAGAGGCTCTG
CCAACCGAAGCCTAAGGTCCAAACCAAGGCAGGGGCCAAAGCTCCAGCTAAGGCCAGGCTTCAGCTCCA
GCTCAGGCTCCAAAGGTGCTCAGGCCCCAAAGGTGCCAGGCCCTGTGAAGGCCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227076 protein sequence
Red=Cloning site Green=Tags(s)

MAKSKNHTTHNQSRKWHNRNGIKKPRSQRYESLKGVDPKFLRNMRFAKKNKGLKMKQANNAKAVSARAE
AIKALVKPQAIKPKMPKGPPLKRLAFIAHPKLGKIRISYMAKGQRLCQPKPKVQTKAGAKAPAKAQASAP
AQAPKGAQAPKGAQAPVKAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:


ACCN: NM_009082

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009082.3](#), [NP_033108.1](#)

RefSeq Size: 710 bp


RefSeq ORF: 483 bp

Locus ID: 19944

UniProt ID: [P47915](#)

Cytogenetics: 9 F1
 MW: 17.6 kDa
 Gene Summary: Component of the large ribosomal subunit.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227076