

Product datasheet for **MR227119**

Nr3c1 (NM_008173) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nr3c1 (NM_008173) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nr3c1
Synonyms:	GR; Grl-1; Grl1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Protein Sequence: >MR227119 representing NM_008173
Red=Cloning site Green=Tags(s)


MDSKESLAPPGRDEVPSLLGRGRGSVMDLYKTLRGGATVKVSASSPSVAAASQADSKQQRILLDFSKGS
ASNAQQQQQQQQQQQQQQQPDL SKAVSL SMGL YMGETETKVMGNDLGYPQQQLGLSSGETDFRLL
EESIANLNRSTSRPENPKSSTPAAGCATPTEKEFPQTHSDPSSEQQNRKSQPGTNGGSKLYTTDQSTFD
ILQDLEFSAGSPGKETNESPWRSDDLIDENLLSPLAGEDDPFLLEGDVNEDCKPLILPDTKPKIQDTGDT
ILSSPSSVALPQVKTEKDDFIELCTPGVIKQEKLGPVYQASFSGTNII GNKMSAISVHGVSTSGGQMYH
YDMNTASLSQQQDQKPVFVNIPIPVGSENWNRCSQSGEDNLTSLGAMNFAGRSVFSNGYSSPGMRPDVS
SPPSSSSTATGPPPKLCLVCSDEASGCHYGVLTCGSCKVFFKRAVEGQHNYLCAGRNDICIIDKIRRNKCP
ACRYRKCLQAGMNL EARKTKKIKGIQATAGVSQDTSENANKTIVPAALPQLTPTLVSLLEVIEPEVLY
AGYDSSVPDSAWRIMTTLNMLGGRQVIAAVKWAKAIPGFRNLHLDDQMTLLQYSWMFLMAFALGWRSYRQ
ASGNLLCFAPDLIINEQRMTLPCMYDQCKHMLFISTELQRLQVSYEEYLCMKTLLLLSSVPKEGLKSQEL
FDEIRMTYIKELGKAIVKREGNSSQNWRFYQLTKLLDSMHDVVENLLSYCFQTFLDKSMSEIEFPEMLAE
IITNQIPKYSNGNIKLLLFHQK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja2271_g03.zip

Restriction Sites: Sgfl-Mlul

Cloning Scheme:



ACCN: NM_008173


ORF Size: 2379 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

- OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
- Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
- RefSeq:** [NM_008173.4](#)
- RefSeq Size:** 6345 bp
- RefSeq ORF:** 2379 bp
- Locus ID:** 14815
- Cytogenetics:** 18 21.09 cM
- MW:** 87.2 kDa

Product images:



Circular map for MR227119