

Product datasheet for **MR227119A1V**

Mouse Nr3c1 (NM_008173) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Nr3c1 (NM_008173) AAV Particle
Tag:	Myc-DDK
Symbol:	Nr3c1
Synonyms:	GR; Grl-1; Grl1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

ORF Nucleotide Sequence:

```
>MR227119 representing NM_008173
Red=Cloning site Blue=ORF Green=Tags(s)
```

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR227119 representing NM_008173 Red=Cloning site Green=Tags(s) MDSKESLAPPGRDEVPSLLGRGRGSVMDLYKTLRGGATVKVSASSPSVAAASQADSKQQRILLDFSKGS ASNAQQQQQQQQQQQQQQQPDLSKAVSLSMGLYMGETETKVMGNDLGYPQQQGLSSGETDFRLL EESIANLNRSTSRPENPKSSTPAAGCATPTEKEFPQTHSDPSSEQQNRKSQPGTNGGSKLYTTDQSTFD ILQDLEFSAGSPGKETNESPWRSDDLIDENLLSPLAGEDDPFLLEGDVNEDCKPLILPDTKPKIQDTGDT ILSSPSSVALPQVKTEKDDFIELCTPGVIKQEKLGPVYCQASFSGTNIIGNKMSAISVHGVSTSGGQMYH YDMNTASLSQQQDQKPVFNVIPPIPVGSENWNRQCQSGEDNLTSLGAMNFAGRSVFSNGYSSPGMRPDVS SPPSSSSTATGPPPKLCLVCSDEASGCHYGVLTCGSCKVFFKRAVEGQHNYLCAGRNDICIIDKIRKNCPC ACRYRKCLQAGMNLEARKTKKKIKGIQQATAGVSQDTSENANKTIVPAALPQLTPTLVSLLEVIEPEVLY AGYDSSVPDSAWRIMTTLNMLGGRQVIAAVKWAKAIPGFRNLHDDQMTLLQYSWMFLMAFALGWRSYRQ ASGNLLCFAPDLIINEQRMTLPCMYDQCKHMLFISTELQRLQVSYEEYLCMKTLSSVPKEGLKSQEL FDEIRMTYIKELGKAIVKREGNSSQNWRFYQLTKLLDSMHDVVENLLSYCFQTFLDKSMSIEFPEMLAE IITNQIPKYSNGNIKKLLFHQK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_008173
ORF Size:	2376 bp
Buffer:	PBS with 0.001% Pluronic F4234
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_008173.3</u>
RefSeq Size:	6345 bp
RefSeq ORF:	2379 bp
Locus ID:	14815
Cytogenetics:	18 21.09 cM
MW:	87.2 kDa