

## Product datasheet for MR227125

### Cd40 (NM\_170704) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cd40 (NM_170704) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cd40
Synonyms:	AI326936; Bp50; GP39; HIGM1; IGM; IMD3; p50; T-BAM; Tnfrsf5; TRAP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR227125 representing NM_170704 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTGTCTTTGCCTCGGCTGTGCGGCTATGGGGCTGCTTGTGACAGCGGTCCATCTAGGGCAGTG  
TTACGTGCAGTGACAAACAGTACCTCCACGATGGCCAGTGTGTGATTTGTCCAGCCAGGAAGCCGACT  
GACAAGCCACTGCACAGCTCTTGAGAAGACCAATGCCACCCATGTGACTCAGGCGAATTCAGCCAG  
TGGAACAGGGAGATTCGCTGTCCACGACAGACTGTGAACCAATCAAGGGCTTCGGGTTAAGAAGG  
AGGGCACCGCAGAATCAGACACTGTCTGTACCTGTAAGGAAGGACAACACTGCACCAGCAAGGATTGCGA  
GGCATGTGCTCAGCACACGCCCTGTATCCCTGGCTTTGGAGTTATGGAGATGGCCACTGAGACCACTGAT  
ACCGTCTGTATCCCTGCCAGTCGGCTTCTTCTCCAATCAGTCATCACTTTTCGAAAAGTGTATCCCT  
GGACAAGCTGTGAGGATAAGAATTGGAGGTCCTACAGAAAGGAACGAGTCAGACTAATGTCATCTGTGG  
TTAAAGTCCCGGATGCGAGCCCTGCTGGTCATTCCTGTCGTGATGGGCATCCTCATCACCATTTTCGGG  
GTGTTTCTATATCAGTGGTCAAGAAACCAAGGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


**Protein Sequence:** >MR227125 representing NM\_170704  
Red=Cloning site Green=Tags(s)

MVSLPRLCALWGCLLTAVHLGQCVTCSDKQYLHDGQCCDLCPGSRLTSHCTALEKTQCHPCDSGEFSAQ  
 WNREIRCHQHRHCEPNQGLRVKKEGTAESDTVCTCKEQHCTSKDCEACAQHTPCIPGFGVMEMATETTD  
 TVCHPCPVGFFSNQSSLFEKCPWTSCEDKNLEVLQKGTSTQTNVICGLKSRMRALLVIPVVMGILITIFG  
 VFLYISGQETKG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_170704

**ORF Size:** 666 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_170704.2</a> , <a href="#">NP_733805.1</a>
<b>RefSeq Size:</b>	1678 bp
<b>RefSeq ORF:</b>	669 bp
<b>Locus ID:</b>	21939
<b>UniProt ID:</b>	<a href="#">P27512</a>
<b>Cytogenetics:</b>	2 85.38 cM
<b>MW:</b>	24.9 kDa
<b>Gene Summary:</b>	Receptor for TNFSF5/CD40LG. Transduces TRAF6- and MAP3K8-mediated signals that activate ERK in macrophages and B cells, leading to induction of immunoglobulin secretion. [UniProtKB/Swiss-Prot Function]