

Product datasheet for **MR227144L4V**

Esr2 (NM_207707) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Esr2 (NM_207707) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Esr2
Synonyms:	ERbeta; ER[b]; Estrb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207707
ORF Size:	1701 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227144).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207707.1 , NP_997590.1
RefSeq Size:	3362 bp
RefSeq ORF:	1704 bp
Locus ID:	13983
UniProt ID:	O08537
Cytogenetics:	12 33.52 cM



[View online »](#)

Gene Summary:

Nuclear hormone receptor. Binds estrogens with an affinity similar to that of ESR1 (ER-alpha), and activates expression of reporter genes containing estrogen response elements (ERE) in an estrogen-dependent manner. May play a role in ovarian follicular growth and maturation. [UniProtKB/Swiss-Prot Function]