

## Product datasheet for **MR227146L3V**

### **Dnajb11 (NM\_026400) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Dnajb11 (NM_026400) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Dnajb11
Synonyms:	1810031F23Rik; ABBP-2; AL024055; Dj9; ERdj3; ERj3p
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_026400
ORF Size:	1074 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227146).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_026400.5</a> , <a href="#">NP_080676.3</a>
RefSeq Size:	2534 bp
RefSeq ORF:	1077 bp
Locus ID:	67838
UniProt ID:	<a href="#">Q99KV1</a>
Cytogenetics:	16 B1



[View online »](#)

**Gene Summary:**

As a co-chaperone for HSPA5 it is required for proper folding, trafficking or degradation of proteins. Binds directly to both unfolded proteins that are substrates for ERAD and nascent unfolded peptide chains, but dissociates from the HSPA5-unfolded protein complex before folding is completed. May help recruiting HSPA5 and other chaperones to the substrate. Stimulates HSPA5 ATPase activity. It is necessary for maturation and correct trafficking of PKD1.[UniProtKB/Swiss-Prot Function]