

Product datasheet for **MR227153A1V**

Mouse Ccl5 (NM_013653) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ccl5 (NM_013653) AAV Particle
Tag: Myc-DDK
Symbol: Ccl5
Synonyms: MuRantes; RANTES; Scya5; SISd; TCP228
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR227153 representing NM_013653
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAAGATCTCTGCAGCTGCCCTCACCATCATCCTCACTGCAGCCGCCCTCTGCACCCCGCACCTGCCT
CACCATATGGCTCGGACACCACTCCCTGCTGCTTTGCCTACCTCTCCCTCGCGCTGCCTCGTGCCACGT
CAAGGAGTATTCTACACCAGCAGCAAGTGCTCCAATCTGCAGTCGTGTTTGTCACTCGAAGGAACCGC
CAAGTGTGTGCCAACCCAGAGAAGAAGTGGGTTCAAGAATACATCAACTATTTGGAGATGAGC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227153 representing NM_013653
Red=Cloning site **Green**=Tags(s)

MKISAAALTIIILTAAALCTPAPASPYGSDTTPCCFAYLSLALPRAHVKEYFYTSSKCSNLAVVFVTRNR
QVCANPEKKWVQEYINYLEMS

TRTRPLE**QKLISEEDLA**AND**ILDYKDDDDKV**

Species: Mouse
Serotype: AAV-2
ACCN: NM_013653
ORF Size: 273 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013653.3</u> , <u>NP_038681.2</u>
RefSeq Size:	579 bp
RefSeq ORF:	276 bp
Locus ID:	20304
UniProt ID:	<u>P30882</u>
Cytogenetics:	11 50.66 cM
MW:	10.5 kDa