

## Product datasheet for **MR227162L4V**

### Timp1 (NM\_011593) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Timp1 (NM_011593) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Timp1
Synonyms:	Clgj; EPA; Timp; TIMP-1; TPA-S1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011593
ORF Size:	615 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227162).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_011593.2</a> , <a href="#">NP_035723.2</a>
RefSeq Size:	902 bp
RefSeq ORF:	618 bp
Locus ID:	21857
UniProt ID:	<a href="#">P12032</a>
Cytogenetics:	X 16.38 cM



[View online »](#)

**Gene Summary:**

Metalloproteinase inhibitor that functions by forming one to one complexes with target metalloproteinases, such as collagenases, and irreversibly inactivates them by binding to their catalytic zinc cofactor. Acts on MMP1, MMP2, MMP3, MMP7, MMP8, MMP9, MMP10, MMP11, MMP12, MMP13 and MMP16. Does not act on MMP14 (By similarity). Also functions as a growth factor that regulates cell differentiation, migration and cell death and activates cellular signaling cascades via CD63 and ITGB1. Plays a role in integrin signaling.  
[UniProtKB/Swiss-Prot Function]