

Product datasheet for **MR227167L3V**

Irs2 (NM_001081212) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Irs2 (NM_001081212) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Irs2
Synonyms:	Irs-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001081212
ORF Size:	3963 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227167).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001081212.1
RefSeq Size:	4015 bp
RefSeq ORF:	3966 bp
Locus ID:	384783
UniProt ID:	P81122
Cytogenetics:	8 5.35 cM
Gene Summary:	May mediate the control of various cellular processes by insulin.[UniProtKB/Swiss-Prot Function]



[View online »](#)