

Product datasheet for **MR227168L3V**

Pax7 (NM_011039) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pax7 (NM_011039) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pax7
Synonyms:	Pax-7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011039
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227168).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011039.2 , NP_035169.1
RefSeq Size:	4275 bp
RefSeq ORF:	1512 bp
Locus ID:	18509
UniProt ID:	P47239
Cytogenetics:	4 70.83 cM
Gene Summary:	Transcription factor playing a role in myogenesis through regulation of muscle precursor cells proliferation.[UniProtKB/Swiss-Prot Function]



[View online »](#)