

Product datasheet for **MR227178**

Irak1 (NM_001177974) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Irak1 (NM_001177974) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Irak1
Synonyms:	AA408924; Il1rak; IRAK; IRAK-1; IRAK1-S; IRAK1b; mPLK; Plpk
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR227178 representing NM_001177974
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGGGGGGCCGGGCCCGGGGAGCCTGTGGTTCCAGGCGCCAGCACTTCTTGTACGAGGTGCCAC
CCTGGGTTATGTGCCGTTCTACAAAGTGATGGACGCCCTGGAGCCCGCGACTGGTGCCAGTTCGCGGC
CTTGATCGTGCAGCAGCAGAGCTGCGGCTGTGCGAGCGCTCCGAGCAGCGCACAGCCAGTGTCTG
TGGCCCTGGATCAACCGCAACGCGCGCGTAGTGACCTCGTTCACATCCTCACGCACCTGCAGCTGCTGC
GTGCGAGGGACATCATCACAGCCTGGCACCCCTCTGCCCCGTTGTGCCCAAGCACCGCTGCCCAAG
GCCAGCAGCATCTGCAAGGCTCTGAGGCCGGGACTGGAGCCCCGGAAATTGCAGTCTCTGCCTCC
ACCTTCTCTCCCGAGCTTTCCAGGCTCCAGACCCATTCTGAGTCAGAGCTCTCCAGGTTCCACTCC
CTGTTTCCCTCGGGCCACCCTACCCTTTCAGCCCTTCTCCACCAAGTCAAGCCAGAGAGCCAGT
GTCTGGCTCCAAGAGCCCATCCCTCCCGTTTTGCTGGCCCTTCTGTGAGATTTCCAAGGCACCTGC
AACTTCTCTGAAGAACTCAGATTGGAGAGGGTGGTTTTGGATGTGTGTACCGAGCAGTCATGAGAAATA
CTACATATGCTGTGAAGAGACTGAAGGAGGAAGCTGACCTAGAGTGGACTATGGTGAAACAGAGTCTT
AACAGAGGTGGAACAGCTATCAAGGTTTCGTACCCAAATATCGTAGACTTTGCTGGCTACTGTGCAGAG
AGTGGCTTATACTGCCTTGTATGGCTTCTTGCCCAATGGCTCCTTGGAGGATCAGTCCACCTCAGA
CCCAGGCTGCTCCCACTTTCTGGCCTCAACGACTGGACATTTCTTGGGCACAGCCCGGCTATTCA
GTTTTTACATCAGGATAGCCCCAGCCTTATCCATGGAGACATCAAGAGTTCTAACGTGCTTCTGGATGAG
AGACTGATGCCAAGCTGGGAGACTTTGGCCTGGCTCGTTTTCAGCCGCTTTCGGGGGCCAAAGCAAGC
AGAGCAGTACTGTGGCCGGACTTCCACAGTTCGAGGTACCCTCGCTACTTGCCTGAGGATACATCAA
GACAGGCCGACTGGCTGTGGACACCGATACCTTCAGCTTTGGGGTGGTAATACTGGAGACCCTTGCTGGT
CAGAGGGCTGTGAGGACACAAGGTGCAAAGACCAAGTATTTGAAAGACCTGATTGAAGATGAGGCTGAAG
AGGCTGGAGTGACCCTGAAAAGCACCCAGCCTACTCTGTGGTGGGTGTAGCCACGGATGCTTGGCTGC
TCCAATTGCTGCCAGATCTATAAGAAGCACCTGGACTCCAGACCTGGGCCCTGCCACCCAGTTGGGC
CTGGCCCTGGCTCAACTAGCTTGTGCTGCATGCACCGTCCGGCCAAAGAGGGCCCCCATGACCCAGG
TATACAAGAGACTAGAAGGGCTTCAGGCAGGGCTCCCTGGGAGCTAGAGGTTGCCGGCCATGGCTCCCC
TTCCCCACAGGAGAACTCTACATGTCTACCCTGGCAGTGCCCAGAGTGGGGATGAACCATGGCAGCCT
CTAGTAGTGACCACAAGAGCCCCAGCCAGGCTGCCAGCAACTCCAGAGAAGTCCAACAGCCAGTGG
AAAGTGATGAGAGTGTTCCCGGCCCTCTGCTACCCTGCATTCTGGCACTTGACTCCAGGTTCCACCC
AAGCCCAGCGTCTTCAGAGAGGCTAGCTGTACCCAGGGAGGCACTACCAGAGAATCAAGTGTGAGGAGT
AGCCCAGGCTTCCAGCCTACAACCATGGAAGGCTCACCCACGGGCAGCTCATCCCTGCTGTGCATCAGAAC
CACCACAGATCATCATCAACCCAGCCGACAGAAGATGGTACAAAAGCTGGCTCTTTATGAAGAAGGGGT
CTTGGATAGCCTGCAACTGCTGTGCATCAGGCTTTTCCAGGCTTGGATTTAGAACCTGAAAAGAGCCAG
GGACCTGAAGAAAGTGATGAATTCCAGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227178 representing NM_001177974
 Red=Cloning site Green=Tags(s)


```
MAGGPGPEVPVPGAQHFLYEVPWVMCRFYKVMdalePADWCQFAALIVRDQTELRLCERSEQRTASVL
WPWINRNARVADLVHILTHLQLLRARDIITAWHPPAPVPPSTAAPRPSSISAGSEAGDWSPRKLQSSAS
TFLSPAFPGSQTHSESELLQVPLPVSLGPPLPSSAPSSTKSSPESPVSGLQRAHPSPF CWPFC EISQGTC
NFSEELRIGEGGFGCVYRAVMRNTTYAVKRLKEEADLEWTVKQSFLEVEQLSRFRHPNIVDFAGYCAE
SGLYCLVYGFLPNGSLEDQLHLQTQACSPLSWPQRLDILLGTARAIQFLHQDSPSLIHGDIKSSNVLLDE
RLMPKLGDFLARFSRFAGAKASQSSTVARTSTVRGTLAYLPEEYIKTGR LAVDTDTFSFGVVILETLAG
QRAVRTQGAKT KYLKDLEDEAEAGVTLKSTQPTLWVGVATDAAWAAIAAQIYKHLDSRPGPCPPQLG
LALAQLACCCMHRRAKKRPPMTQVYKRLEGLQAGPPWELEVAGHGSPSPQENSYMSTTGSQAQSGDEPWQP
LVVTTTRAPAQAQQLRSPNQPVESDESVPGLSATLHSHWLTGSHSPSPASFREASCTQGGTTRESSVRS
SPGFQPTTMEGSPGSSLLSSEPPQII INPARQKMVQKLALYEEGVLDLSQLLSSGFFPGLDLEPEKSQ
GPEESDEFQS
```

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001177974

ORF Size: 2130 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177974.1](#), [NP_001171445.1](#)

RefSeq Size: 2574 bp

RefSeq ORF: 2133 bp

Locus ID: 16179


UniProt ID: [Q62406](#)

Cytogenetics: X 37.61 cM

MW: 77.3 kDa

Gene Summary: Serine/threonine-protein kinase that plays a critical role in initiating innate immune response against foreign pathogens. Involved in Toll-like receptor (TLR) and IL-1R signaling pathways. Is rapidly recruited by MYD88 to the receptor-signaling complex upon TLR activation. Association with MYD88 leads to IRAK1 phosphorylation by IRAK4 and subsequent autophosphorylation and kinase activation. Phosphorylates E3 ubiquitin ligases Pellino proteins (PELI1, PELI2 and PELI3) to promote pellino-mediated polyubiquitination of IRAK1. Then, the ubiquitin-binding domain of IKBKG/NEMO binds to polyubiquitinated IRAK1 bringing together the IRAK1-MAP3K7/TAK1-TRAF6 complex and the NEMO-IKKA-IKKB complex. In turn, MAP3K7/TAK1 activates IKKs (CHUK/IKKA and IKBKB/IKKB) leading to NF-kappa-B nuclear translocation and activation. Alternatively, phosphorylates TIRAP to promote its ubiquitination and subsequent degradation. Phosphorylates the interferon regulatory factor 7 (IRF7) to induce its activation and translocation to the nucleus, resulting in transcriptional activation of type I IFN genes, which drive the cell in an antiviral state. When sumoylated, translocates to the nucleus and phosphorylates STAT3 (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227178