

## Product datasheet for **MR227189**

### **Irak1 (NM\_008363) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Irak1 (NM_008363) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Irak1
Synonyms:	AA408924; Il1rak; IRAK; IRAK-1; IRAK1-S; IRAK1b; mPLK; Plpk
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MR227189 representing NM\_008363  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCCGGGGGGCCGGGCCCGGGGAGCCTGTGGTTCCAGGCGCCAGCACTTCTTGTACGAGGTGCCAC  
 CCTGGGTTATGTGCCGCTTCTACAAAGTGATGGACGCCCTGGAGCCCGCGACTGGTGCCAGTTCGCGGC  
 CTTGATCGTGCGGACCAGACAGAGCTGCGGCTGTGCGAGCGCTCCGAGCAGCGCACAGCCAGTGTCTG  
 TGGCCCTGGATCAACCGCAACGCGCGCGTAGTGACCTCGTTCACATCCTCACGCACCTGCAGTGTCTG  
 GTGCGAGGGACATCATCACAGCCTGGCACCCCTCTGCCCCGTTGTGCCCAAGCACCGCTGCCCAAG  
 GCCCAGCAGCATCTGCAAGGCTCTGAGGCCGGGACTGGAGCCCCGAAATTGCAGTCTCTGCCTCC  
 ACCTTCTCTCCCGAGCTTTCCAGGCTCCAGACCCATTCTGAGTCAGAGCTCTCCAGGTTCCACTCC  
 CTGTTTCCCTCGGGCCACCCTACCCTTTCAGCCCTTCTCCACCAAGTCAAGCCAGAGAGCCAGT  
 GTCTGGCTCCAAGAGCCCATCCCTCCCGTTTTGCTGGCCCTTCTGTGAGATTTCCAAGGCACCTGC  
 AACTTCTCTGAAGAACTCAGATTGGAGAGGGTGGTTTTGGATGTGTGTACCGAGCAGTCATGAGAAATA  
 CTACATATGCTGTGAAGAGACTGAAGGAGGAAGCTGACCTAGAGTGGACTATGGTGAAACAGAGTCTT  
 AACAGAGGTGGAACAGCTATCAAGGTTTCGTACCCAAATATCGTAGACTTTGCTGGCTACTGTGCAGAG  
 AGTGGCTTATACTGCCTTGTATGGCTTCTTGCCCAATGGCTCCTTGGAGGATCAGTCCACCTTCCAGA  
 CCCAGGCTGCTCCCACTTTCTGGCCTCAACGACTGGACATTTCTTGGGCACAGCCCGGCTATTCA  
 GTTTTTACATCAGGATAGCCCCAGCCTTATCCATGGAGACATCAAGAGTTCTAACGTGCTTCTGGATGAG  
 AGACTGATGCCAAGCTGGGAGACTTTGGCCTGGCTCGTTTCAGCCGCTTTCGGGGGCCAAAGCAAGC  
 AGAGCAGTACTGTGGCCGGACTTCCACAGTTCGAGGTACCCTCGCCTACTTGCCTGAGGATACATCAA  
 GACAGGCCGACTGGCTGTGGACACCGATACCTTCAGCTTTGGGGTGGTAATACTGGAGACCCTTGCTGGT  
 CAGAGGGCTGTGAGGACACAAGGTGCAAAGACCAAGTATTTGAAAGACCTGATTGAAGATGAGGCTGAAG  
 AGGCTGGAGTGACCCTGAAAAGCACCCAGCCTACTCTGTGGTGGGTGTAGCCACGGATGCTTGGCTGC  
 TCCAATTGCTGCCAGATCTATAAGAAGCACCTGGACTCCAGACCTGGGCCCTGCCACCCAGTTGGGC  
 CTGGCCCTGGCTCAACTAGCTTGTGCTGCATGCACCGTCCGGCCAAAGAGGGCCCCCATGACCCAGG  
 TATACAAGAGACTAGAAGGGCTTCAGGCAGGGCTCCCTGGGAGCTAGAGGTTGCCGGCCATGGCTCCCC  
 TTCCCCACAGGAGAACTCTACATGTCTACCCTGGCAGTCCCAGAGTGGGGATGAACCATGGCAGCCT  
 CTAGTAGTGACCACAAGAGCCCCAGCCAGGCTGCCAGCAACTCCAGAGAAGTCCAACAGCCAGTGG  
 AAAGTGATGAGAGTGTTCCCGGCCCTCTCTGCTACCCTGCATTCTGGCACTTGACTCCAGGTTCCACCC  
 AAGCCAGCGTCTTTCAGAGAGGCTAGCTGTACCCAGGGAGGCACTACCAGAGAATCAAGTGTGAGGAGT  
 AGCCAGGCTTCCAGCCTACAACCATGGAAGGCTCACCCACGGGCAGCTCATCCCTGCTGTGCATCAGAAC  
 CACCACAGATCATCATCAACCCAGCCGACAGAAGATGGTACAAAAGCTGGCTCTTTATGAAGAAGGGGT  
 CTTGGATAGCCTGCAACTGCTGTGCATCAGGCTTTTCCAGGCTTGGATTTAGAACCTGAAAAGAGCCAG  
 GGACCTGAAGAAAGTGATGAATTCCAGAGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR227189 representing NM\_008363  
 Red=Cloning site Green=Tags(s)


MAGGPGPEVPVPGAQHFLYEVPWVMCRFYKVMdalePADWCQFAALIVRDQTELRlcERSEQRTASVL  
 WPWINRNARVADLVHILTHLQLLRARDIITAWHPPAPVVPSTAAAPRPSSISAGSEAGDWSPRKLQSSAS  
 TFLSPAFFPGSQTHSESELLQVPLPVSLGPPLPSSAPSSTKSSPESPVSGLQRAHPSPF CWPFC EISQGTC  
 NFSEELRIGEGGFGCVYRAVMRNTTYAVKRLKEEADLEWTMVKQSFLTEVEQLSRFRHPNIVDFAGYCAE  
 SGLYL VYGFLPNGSLEDQLHLQTQACSPLSWPQRLDILLGTARAIQFLHQDSPSLIHGDIKSSNVLLDE  
 RLMPKLGDFLARFSRFAGAKASQSSTVARTSTVRGTLAYLPEEYIKTGR LAVDTDTFSFGVVILETLAG  
 QRAVRTQGAKTKYLKDLIEDEAEAEAGVTLKSTQPTLWVG VATDAWAAPIAAQIYKHLDSRPGPCPPQLG  
 LALAQLACCCMHRRAKKRPMTQVYKRLEGLQAGPPWELEVAGHGSPSPQENSYMSTTGSQAQSGDEPWQP  
 LVVTTTRAPAAQQLQRSPNQPVESDESVPGLSATLHSHWLTGSHSPSPASFREASCTQGTTRESSVRS  
 SPGFQPTTMEGSP TGS SLLSSEPPQII INPARQKMQK LALYEEGVLDSLQLLSSGFFPGLDLEPEKSQ  
 GPEESDEFQS

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2166\\_d04.zip](https://cdn.origene.com/chromatograms/ja2166_d04.zip)

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_008363

**ORF Size:** 2130 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_008363.2](#), [NP\\_032389.2](#)

**RefSeq Size:** 3831 bp

**RefSeq ORF:** 2133 bp

**Locus ID:** 16179


**UniProt ID:** [Q62406](#)

**Cytogenetics:** X 37.61 cM

**MW:** 77.7 kDa

**Gene Summary:**

Serine/threonine-protein kinase that plays a critical role in initiating innate immune response against foreign pathogens. Involved in Toll-like receptor (TLR) and IL-1R signaling pathways. Is rapidly recruited by MYD88 to the receptor-signaling complex upon TLR activation. Association with MYD88 leads to IRAK1 phosphorylation by IRAK4 and subsequent autophosphorylation and kinase activation. Phosphorylates E3 ubiquitin ligases Pellino proteins (PELI1, PELI2 and PELI3) to promote pellino-mediated polyubiquitination of IRAK1. Then, the ubiquitin-binding domain of IKBKG/NEMO binds to polyubiquitinated IRAK1 bringing together the IRAK1-MAP3K7/TAK1-TRAF6 complex and the NEMO-IKKA-IKKB complex. In turn, MAP3K7/TAK1 activates IKKs (CHUK/IKKA and IKBKB/IKKB) leading to NF-kappa-B nuclear translocation and activation. Alternatively, phosphorylates TIRAP to promote its ubiquitination and subsequent degradation. Phosphorylates the interferon regulatory factor 7 (IRF7) to induce its activation and translocation to the nucleus, resulting in transcriptional activation of type I IFN genes, which drive the cell in an antiviral state. When sumoylated, translocates to the nucleus and phosphorylates STAT3 (By similarity). [UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for MR227189