

## Product datasheet for **MR227196L3V**

### Gfi1 (NM\_010278) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gfi1 (NM_010278) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gfi1
Synonyms:	AW495828; Gfi-1; Pal-1; Pal1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010278
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_010278.2</a> , <a href="#">NP_034408.1</a>
RefSeq Size:	2862 bp
RefSeq ORF:	1272 bp
Locus ID:	14581
UniProt ID:	<a href="#">P70338</a>
Cytogenetics:	5 52.23 cM



[View online »](#)

**Gene Summary:**

Transcription repressor essential for hematopoiesis. Functions in a cell-context and development-specific manner. Binds to 5'-TAAATCAC[AT]GCA-3' in the promoter region of a large number of genes. Component of several complexes, including the EHMT2-GFI1-HDAC1, AJUBA-GFI1-HDAC1 and RCOR-GFI-KDM1A-HDAC complexes, that suppress, via histone deacetylase (HDAC) recruitment, a number of genes implicated in multilineage blood cell development. Regulates neutrophil differentiation, promotes proliferation of lymphoid cells, and is required for granulocyte development. Mediates, together with U2AF1L4, the alternative splicing of CD45 and controls T-cell receptor signaling. Regulates the endotoxin-mediated Toll-like receptor (TLR) inflammatory response by antagonizing RELA. Cooperates with CBFA2T2 to regulate ITGB1-dependent neurite growth. Controls cell-cycle progression by repressing CDKN1A/p21 transcription in response to TGFB1 via recruitment of GFI1 by ZBTB17 to the CDKN1A/p21 and CDKN1B promoters. Required for the maintenance of inner ear hair cells.[UniProtKB/Swiss-Prot Function]