

Product datasheet for MR227241L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Abcg1 (NM_009593) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Abcg1 (NM_009593) Mouse Tagged ORF Clone Lentiviral Particle

Symbol: Abcg⁻

Synonyms: Abc; Abc8; AW413978; Whit; White

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_009593 **ORF Size:** 1998 bp

ORF Nucleotide

- - - - -

Sequence:

The ORF insert of this clone is exactly the same as(MR227241).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 009593.2, NP 033723.1

RefSeq Size: 5821 bp
RefSeq ORF: 2001 bp
Locus ID: 11307
UniProt ID: Q64343

Cytogenetics: 17 A3.3







Gene Summary:

The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. The human protein is involved in macrophage cholesterol efflux and may regulate cellular lipid homeostasis in other cell types. [provided by RefSeq, Jul 2008]