

Product datasheet for **MR227242**

Hoxb9 (NM_008270) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hoxb9 (NM_008270) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Hoxb9
Synonyms: Hox-2.5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR227242 representing NM_008270
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCATTTCTGGGACGCTTAGCAGCTATTATGTCGACTCGATCATAAGTCACGAGAGCGAGGACGCGC
CTCCAGCCAAGTTTCCTTCGGGCCAGTACGCGAACCCGCGGCAGCCGGGCCACGCCGAGCACCTGGACTT
CCCCTCGTGCAGCTTCCAGCCCAAAGCGCCGTGTTTCGGCGCCTCCTGGGCGCCGCTGAGCCCGCACGCC
TCGGGGAGCCTGCCGTCTGTCTACCACCCTACCTGCAGCCCCAGGGCGCGCCGGCGCCGAAAGCAGGT
ACCTCCGCACCTGGCTAGAGCCGGCGCCGCGCGAGGCGCCCGGGGAGGGCCAGGCGCCGCTGAA
GGCGGAGCCGCTGCTGGGCGCGCCGGGGAGCTGCTCAAACAGGGCACGCCCGAGTACAGTTTGAAACT
TCGGCGGGCAGGGAGGCTGTCCTGTCTAATCAAAGAGCTGGCTACGGGGACAATAAAATTTGCGAAGGAA
GCGAGGACAAAGAGAGGCCGGATCAAACCAACCCCTCTGCCAACTGGCTGCACGCCCGCTCTTCCCGGAA
AAAGCGCTGTCCCTACACCAAATACCAGACGCTGGAGCTAGAGAAGGAGTTTCTGTTCAATATGTACCTC
ACCAGGGACCGCAGGCATGAAGTGGCCAGACTCCTCAATCTGAGTGAGAGACAAGTCAAATCTGGTTTC
AGAACCGCGGATGAAAATGAAGAAGATGAATAAGGAGCAGGGCAAAGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008270.2](#), [NP_032296.2](#)

RefSeq Size: 2711 bp

RefSeq ORF: 753 bp

Locus ID: 15417


UniProt ID: [P20615](#)

Cytogenetics: 11 59.81 cM

MW: 28.5 kDa

Gene Summary: Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227242