

Product datasheet for **MR227247L3V**

Lat2 (NM_020044) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lat2 (NM_020044) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Lat2
Synonyms:	AW125574; LAB; NTAL; Wbscr5; Wbscr15; WSCR5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020044
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020044.2 , NP_064428.1
RefSeq Size:	1523 bp
RefSeq ORF:	612 bp
Locus ID:	56743
UniProt ID:	Q9JHL0
Cytogenetics:	5 G2



[View online »](#)

Gene Summary:

Involved in FCER1 (high affinity immunoglobulin epsilon receptor)-mediated signaling in mast cells. May also be involved in BCR (B-cell antigen receptor)-mediated signaling in B-cells and FCGR1 (high affinity immunoglobulin gamma Fc receptor I)-mediated signaling in myeloid cells. Couples activation of these receptors and their associated kinases with distal intracellular events through the recruitment of GRB2.[UniProtKB/Swiss-Prot Function]