

Product datasheet for MR227254

Mfap2 (NM_008546) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Mfap2 (NM_008546) Mouse Tagged ORF Clone

Tag: Myc-DDK

Symbol: Mfap2

Synonyms: AI893631; Magp; Magp1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >MR227254 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

**ATGAGAGCCGCTGCCTCTTCTGCTGTTCATGCCTGGCCTGCTGGCTCAGGGCCAATATGACCTGGATC
CTCTCCCCCATCCCGACCATGTCCAGTACAACCCTATGGCGACCAGATTGACAACGCAGACTACTA
TGACTACCAAGAAGTGAGTCCTCGGACCCCTGAAGAGCAGTCCAGTCCCAGCAGCAAGTTCAACAGGAA
GTCATCCAGCCCTACCCAGAACCAGCAGCTGCAGGGGACCTGGAGACTGAGCCTACCGAGCCTGGCC
CTTTGACTGCCGGAAGAACAGTACCCGTGTACCCGCCTCTACTCCATCCACAAGCCTTGCAAACAGTG
TCTCAACGAGGTCTGCTTCTACAGCCTCCGCCGAGTGTATGTGGTCAACAAGGAAATATGTGTCCGCACT
GTCTGTGCCACGAGGAACTTCTCCGAGCTGACCTGTGCCGAGACAAGTTCTCCAAGTGTGGCGTATGG
CCGTCACTGGCCTGTGCCAATCTGTGGCTGCCTCTGTGCCAGGAGCTGTGGGGCTGC**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence: >MR227254 protein sequence
Red=Cloning site Green=Tags(s)

MRAACLFLFMPGLLAQGQYDLDPFPFDPHVQYNHYGDQIDNADYYDYQEVSPRTPPEEQFSQQVQQE
VIPAPTPAAPAAAGDLETEPTEPGLDPCREEQYPCRTLYSIHKPKQCLNEVCFYSLRRVYVVNKEICVRT
VCAHEELLRADLCRDKFSKCGVMAVSGLCQSVAAASCARSCGGC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



View online »

Cloning Scheme:


ACCN: NM_008546

ORF Size: 552 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_008546.2](#)

RefSeq Size: 1084 bp

RefSeq ORF: 558 bp

Locus ID: 17150


UniProt ID: [P55002](#)

Cytogenetics: 4 D3

MW: 20.6 kDa

Gene Summary: Component of the elastin-associated microfibrils.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR227254