

## Product datasheet for MR227258

### Hsd11b1 (NM\_008288) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hsd11b1 (NM_008288) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hsd11b1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR227258 representing NM_008288 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCAGTTATGAAAAATTACCTCCTCCCGATCCTGGTGCTCTTCTGGCCTACTACTACTATTCTACAA  
ATGAAGAGTTCAGACCAGAAATGCTCCAGGAAAGAAAGTGATTGTCCTGGGCCAGCAAAGGGATTGG  
AAGAGAAATGGCATATCATCTGTCAAAAATGGGAGCCCATGTGGTATTGACTGCCAGGTCGGAGGAAGGT  
CTCCAGAAGGTAGTGTCTCGCTGCCTTGAACCTCGGAGCAGCCTCTGCTCACTACATTGCTGGCACTATGG  
AAGACATGACATTTGCGGAGCAATTTATTGTCAAGGCGGAAAGCTCATGGGCGGACTGGACATGCTTAT  
TCTAAACCACATCACTCAGACCTCGCTGTCTCTTCCATGACGACATCCACTCTGTGCGAAGAGTCATG  
GAGGTCAACTTCCCTCAGCTACGTGGTCATGAGCACAGCCGCCTTGCCCATGCTGAAGCAGAGCAATGGCA  
GCATTGCCGTCATCTCCTCCTGGCTGGGAAAATGACCCAGCCTATGATTGCTCCCTACTCTGCAAGCAA  
GTTTGCTCTGGATGGTTCTTTCCACCATTAGAACAGAACTCTACATAACCAAGGTCAACGTGTCCATC  
ACTCTCTGTGCTTGGCCTCATAGACACAGAAACAGCTATGAAGGAAATCTCTGGGATAATTAACGCC  
AAGCTTCTCCAAGGAGGAGTGCGCCCTGGAGATCATCAAAGGCACAGCTCTACGCAAAGCGAGGTGTA  
CTATGACAAATCGCCTTTGACTCCAATCCTGCTGGGAACCCAGGAAGGAAGATCATGGAATTTTTTCA  
TTACGATATTATAATAAGGACATGTTTGAAGTAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR227258 representing NM\_008288  
Red=Cloning site Green=Tags(s)

MAVMKNYLLPILVFLAYYYYSTNEEFRPEMLQGKKVIVTGASKGIGREMAHYLSKMGAHVVL TARSEEG  
 LQKVVSRCLELGAASAHYIAGTMEDMTFAEQFIVKAGKLMGGLDMLILNHITQTSLSLFHDDIHSVRRVM  
 EVNFLSYVVMSTAALPMLKQSNQSI AVISSLAGKMTQPMIAPYSASKFALDGGFFSTIRTELYITKVNYSI  
 TLCVLGLIDTETAMKEISGIINAQASPKEECALEIIKGTALRKSEVYDYKSP LTPILLGNPGRKIMEFFS  
 LRYYNKDMFVSN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mm9023\\_g09.zip](https://cdn.origene.com/chromatograms/mm9023_g09.zip)

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_008288


**ORF Size:** 876 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_008288.2</a></u> , <u><a href="#">NP_032314.2</a></u>
<b>RefSeq Size:</b>	1381 bp
<b>RefSeq ORF:</b>	879 bp
<b>Locus ID:</b>	15483
<b>UniProt ID:</b>	<u><a href="#">P50172</a></u>
<b>Cytogenetics:</b>	1 H6
<b>MW:</b>	32.8 kDa
<b>Gene Summary:</b>	Catalyzes reversibly the conversion of cortisol to the inactive metabolite cortisone. [UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR227258