

Product datasheet for MR227281L2V

II6 (NM_031168) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles Product Name: II6 (NM 031168) Mouse Tagged ORF Clone Lentiviral Particle Symbol: 116 II: II-6 Synonyms: **Mammalian Cell** None Selection: Vector: pLenti-C-mGFP (PS100071) mGFP Tag: ACCN: NM 031168 ORF Size: 633 bp The ORF insert of this clone is exactly the same as(MR227281). **ORF** Nucleotide Sequence: **OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info **OTI** Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. **RefSeq:** NM 031168.1 **RefSeq Size:** 1087 bp **RefSeq ORF:** 636 bp Locus ID: 16193 **UniProt ID:** P08505 Cytogenetics: 5 15.7 cM



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



Gene Summary:This gene encodes a member of the interleukin family of cytokines that have important
functions in immune response, hematopoiesis, inflammation and the acute phase response.
The ectopic overexpression of the encoded protein in mice results in excessive plasma cells in
circulation, leading to death. Mice lacking the encoded protein exhibit abnormalities in
hepatic acute phase response, some immune mechanisms, bone resorption in response to
estrogen, liver regeneration and wound healing. Alternative splicing results in multiple
transcript variants encoding different isoforms. [provided by RefSeq, Sep 2015]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US