

## Product datasheet for **MR227284L1V**

### **Cdkn2a (NM\_001040654) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cdkn2a (NM_001040654) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cdkn2a
Synonyms:	Arf; ARF-INK4a; INK4a-ARF; Ink4a/Arf; MTS1; p16; p16(INK4a); p16INK4a; p19p19ARF; Pctr1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001040654
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227284).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001040654.1</a> , <a href="#">NP_001035744.1</a>
RefSeq Size:	850 bp
RefSeq ORF:	507 bp
Locus ID:	12578
UniProt ID:	<a href="#">P51480</a>
Cytogenetics:	4 42.15 cM
Gene Summary:	Acts as a negative regulator of the proliferation of normal cells by interacting strongly with CDK4 and CDK6. This inhibits their ability to interact with cyclins D and to phosphorylate the retinoblastoma protein.[UniProtKB/Swiss-Prot Function]



[View online »](#)