

## Product datasheet for **MR227295L3V**

### Tbx20 (NM\_020496) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbx20 (NM_020496) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tbx20
Synonyms:	9430010M06Rik; AL022859; Tbx12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020496
ORF Size:	891 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227295).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020496.3</a> , <a href="#">NP_065242.1</a>
RefSeq Size:	1800 bp
RefSeq ORF:	894 bp
Locus ID:	57246
UniProt ID:	<a href="#">Q9ES03</a>
Cytogenetics:	9 A4
Gene Summary:	Acts as a transcriptional activator and repressor required for cardiac development and may have key roles in the maintenance of functional and structural phenotypes in adult heart. [UniProtKB/Swiss-Prot Function]



[View online »](#)