

Product datasheet for **MR227304L4V**

Esr1 (NM_007956) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Esr1 (NM_007956) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Esr1
Synonyms:	E; ER; ER-; ER-alpha; ERa; ERalpha; ER[]; Es; ESR; Estr; Estra; Nr; Nr3a1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_007956
ORF Size:	1797 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227304).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_007956.4
RefSeq Size:	4319 bp
RefSeq ORF:	1800 bp
Locus ID:	13982
UniProt ID:	P19785
Cytogenetics:	10 2.03 cM



[View online »](#)

Gene Summary:

This gene encodes an estrogen receptor, a member of the nuclear hormone family of intracellular receptors. The encoded protein, activated by the sex hormone estrogen, is a transcription factor composed of several domains important for hormone binding, DNA binding, and activation of transcription. Estrogen and its receptors are essential for sexual development and reproductive function, but also play a role in other tissues such as bone. Similar genes in human have been implicated in pathological processes including breast cancer, endometrial cancer, and osteoporosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014]