

Product datasheet for **MR227308L4V**

Flt3l (NM_013520) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Flt3l (NM_013520) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Flt3l
Synonyms:	Flt3lg; Ly72L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013520
ORF Size:	696 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227308).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013520.3 , NP_038548.3
RefSeq Size:	1075 bp
RefSeq ORF:	699 bp
Locus ID:	14256
UniProt ID:	P49772
Cytogenetics:	7 29.14 cM
Gene Summary:	Stimulates the proliferation of early hematopoietic cells by activating FLT3. Synergizes well with a number of other colony stimulating factors and interleukins.[UniProtKB/Swiss-Prot Function]



[View online »](#)