

Product datasheet for **MR227344L1V**

Mst1r (NM_009074) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mst1r (NM_009074) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mst1r
Synonyms:	CD136; CDw136; Fv; Fv-; Fv-2; Fv2; PTK; PTK8; Ron; ST; STK
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_009074
ORF Size:	4134 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227344).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_009074.2 , NP_033100.2
RefSeq Size:	4706 bp
RefSeq ORF:	4137 bp
Locus ID:	19882
UniProt ID:	Q62190
Cytogenetics:	9 F1



[View online »](#)

Gene Summary:

This gene encodes a precursor protein that is proteolytically cleaved to yield an alpha chain and a beta chain which form a membrane-spanning heterodimer. The encoded protein belongs to a family of cell-surface receptor tyrosine kinases involved in signaling from the cell surface to the intracellular environment. The binding of the encoded protein to its ligand, macrophage-stimulating protein, mediates several biological activities including wound healing, tumor immunity, macrophage activation and hematopoiesis as well as cell growth, motility, survival and adhesion. The protein encoded by this gene also functions in early development and the macrophage-mediated inflammatory response. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]