

Product datasheet for **MR227350L4V**

Mmp7 (NM_010810) Mouse Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	Mmp7 (NM_010810) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mmp7
Synonyms:	MA; MAT; MMP-7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010810
ORF Size:	801 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227350).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010810.4 , NP_034940.2
RefSeq Size:	878 bp
RefSeq ORF:	804 bp
Locus ID:	17393
Cytogenetics:	9 2.46 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the matrix metalloproteinase family of extracellular matrix-degrading enzymes that are involved in tissue remodeling, wound repair, progression of atherosclerosis and tumor invasion. The encoded preproprotein undergoes proteolytic processing to generate a mature, zinc-dependent endopeptidase enzyme. Mice lacking the encoded protein are deficient in functional cryptdins. This gene is located in a cluster of other matrix metalloproteinase genes on chromosome 9. [provided by RefSeq, Feb 2016]

Product images:

[MR227350L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with MR227350L4V particle to overexpress human Mmp7-mGFP fusion protein.