

Product datasheet for **MR227365L3V**

Mark3 (NM_022801) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mark3 (NM_022801) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mark3
Synonyms:	1600015G02Rik; A430080F22Rik; C-TAK1; Emk2; ETK-1; ETK1; MPK10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022801
ORF Size:	2187 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227365).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022801.4 , NP_073712.2
RefSeq Size:	3343 bp
RefSeq ORF:	2190 bp
Locus ID:	17169
UniProt ID:	Q03141
Cytogenetics:	12 F1



[View online »](#)

Gene Summary:

Serine/threonine-protein kinase. Involved in the specific phosphorylation of microtubule-associated proteins for MAPT/TAU, MAP2 and MAP4. Phosphorylates CDC25C. Regulates localization and activity of some histone deacetylases by mediating phosphorylation of HDAC7, promoting subsequent interaction between HDAC7 and 14-3-3 and export from the nucleus. Negatively regulates the Hippo signaling pathway and antagonizes the phosphorylation of LATS1. Cooperates with DLG5 to inhibit the kinase activity of STK3/MST2 toward LATS1.[UniProtKB/Swiss-Prot Function]