

## Product datasheet for **MR227367L4V**

### Nefm (NM\_008691) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Nefm (NM_008691) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Nefm
Synonyms:	Nef3; NF-M; NF160; NF165; Nfm
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008691
ORF Size:	2544 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227367).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_008691.2</a> , <a href="#">NP_032717.2</a>
RefSeq Size:	3244 bp
RefSeq ORF:	2547 bp
Locus ID:	18040
UniProt ID:	<a href="#">P08553</a>
Cytogenetics:	14 D1
Gene Summary:	Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber.[UniProtKB/Swiss-Prot Function]



[View online »](#)