

## Product datasheet for **MR227376L3V**

### Rapgef3 (NM\_001177811) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Rapgef3 (NM_001177811) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rapgef3
Synonyms:	2310016P22Rik; 9330170P05Rik; Epac; Epac1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177811
ORF Size:	2727 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227376).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177811.1</a> , <a href="#">NP_001171282.1</a>
RefSeq Size:	3774 bp
RefSeq ORF:	2730 bp
Locus ID:	223864
UniProt ID:	<a href="#">Q8VCC8</a>
Cytogenetics:	15 F1



[View online »](#)

**Gene Summary:**

Guanine nucleotide exchange factor (GEF) for RAP1A and RAP2A small GTPases that is activated by binding cAMP. Through simultaneous binding of PDE3B to RAPGEF3 and PIK3R6 is assembled in a signaling complex in which it activates the PI3K gamma complex and which is involved in angiogenesis. Plays a role in the modulation of the cAMP-induced dynamic control of endothelial barrier function through a pathway that is independent on Rho-mediated signaling. Required for the actin rearrangement at cell-cell junctions, such as stress fibers and junctional actin (By similarity).[UniProtKB/Swiss-Prot Function]