

Product datasheet for **MR227390L3V**

Havcr1 (NM_001166631) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Havcr1 (NM_001166631) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Havcr1
Synonyms:	AI503787; KIM-1; TIM-1; Tim1; Timd1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001166631
ORF Size:	849 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227390).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001166631.1 , NP_001160103.1
RefSeq Size:	2001 bp
RefSeq ORF:	849 bp
Locus ID:	171283
UniProt ID:	Q5QNS5
Cytogenetics:	11 B1.1
Gene Summary:	May play a role in T-helper cell development and the regulation of asthma and allergic diseases. Receptor for TIMD4. May play a role in kidney injury and repair (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)