

Product datasheet for **MR227410L2V**

Tnfrsf9 (NM_001077508) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnfrsf9 (NM_001077508) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfrsf9
Synonyms:	4-1BB; A930040111Rik; AA408498; AI325004; Cd137; CDw137; ILA; Ly63
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001077508
ORF Size:	633 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227410).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001077508.1 , NP_001070976.1
RefSeq Size:	1999 bp
RefSeq ORF:	636 bp
Locus ID:	21942
Cytogenetics:	4 81.52 cM
Gene Summary:	Receptor for TNFSF9/4-1BBL. Possibly active during T cell activation.[UniProtKB/Swiss-Prot Function]



[View online »](#)