

## Product datasheet for **MR227415L4V**

### Irak3 (NM\_028679) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Irak3 (NM_028679) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Irak3
Synonyms:	4833428C18Rik; AI563835; IRAK-M
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_028679
ORF Size:	1827 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227415).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_028679.3</a>
RefSeq Size:	2969 bp
RefSeq ORF:	1830 bp
Locus ID:	73914
UniProt ID:	<a href="#">Q8K4B2</a>
Cytogenetics:	10 D2
Gene Summary:	Inhibits dissociation of IRAK1 and IRAK4 from the Toll-like receptor signaling complex by either inhibiting the phosphorylation of IRAK1 and IRAK4 or stabilizing the receptor complex. [UniProtKB/Swiss-Prot Function]



[View online »](#)