

Product datasheet for MR227434

Mbp (NM_001025259) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mbp (NM_001025259) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mbp
Synonyms: C76307; goll; goll-mbp; Hmb; Hmbpr; jv; jve; mld; R75289; shi
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR227434 representing NM_001025259
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

**ATGGCATCACAGAAGAGACCCTCACAGCGATCCAAGTACCTGGCCACAGCAAGTACCATGGACCATGCCA
 GGCATGGCTTCTCCAAGGCACAGAGACACGGGCATCCTTGACTCCATCGGGCGCTTCTTAGCGGTGA
 CAGGGGTGCGCCAAGCGGGGCTCTGGCAAGGACTCACACACGAGAACTACCCATTATGGCTCCCTGCC
 CAGAAGTCGCAGCACGGCCGGACCAAGATGAAAACCCAGTAGTCCATTTCTTCAAGAACATTGTGACAC
 CTGGAACACCACCTCCATCCAAGGGAAGGGGAGAGGCCTGTCCCTCAGCAGATTTAGCTGGGGAGGAAG
 AGACAGCCGCTCTGGATCTCCATGGCGAGACGC**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR227434 representing NM_001025259
 Red=Cloning site Green=Tags(s)

MASQKRPSQRSKYLATASTMDHARHGFLPRHRDTGILDSIGRFFSGDRGAPKRGSGKDSHTRTTHYGLP
 QKSQHGRTQDENPVVHFFKNIIVTPRTPPPSQGKGRGLSLSRFSWGGDRSRSRSGSPMARR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001025259.2</u> , <u>NP_001020430.1</u>
RefSeq Size:	1985 bp
RefSeq ORF:	387 bp
Locus ID:	17196
UniProt ID:	<u>P04370</u>
Cytogenetics:	18 55.84 cM
MW:	14.7 kDa
Gene Summary:	<p>The protein encoded by the classic Mbp gene is a major constituent of the myelin sheath of oligodendrocytes and Schwann cells in the nervous system. However, Mbp-related transcripts are also present in the bone marrow and the immune system. These mRNAs arise from the long Mbp gene (otherwise called "Golli-Mbp") that contains 3 additional exons located upstream of the classic Mbp exons. Alternative splicing from the Golli and the Mbp transcription start sites gives rise to 2 sets of Mbp-related transcripts and gene products. The Golli mRNAs contain 3 exons unique to Golli-Mbp, spliced in-frame to 1 or more Mbp exons. They encode hybrid proteins that have N-terminal Golli aa sequence linked to Mbp aa sequence. The second family of transcripts contain only Mbp exons and produce the well characterized myelin basic proteins. This complex gene structure is conserved among species suggesting that the Mbp transcription unit is an integral part of the Golli transcription unit and that this arrangement is important for the function and/or regulation of these genes. Mutation of the Mbp gene is associated with the 'shiverer' and 'myelin deficient' phenotypes in mouse. [provided by RefSeq, Jul 2008]</p>