

Product datasheet for **MR227441L3V**

Gemin2 (NM_025656) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gemin2 (NM_025656) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gemin2
Synonyms:	gemin-2; S; Sip1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025656
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227441).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025656.4 , NP_079932.2
RefSeq Size:	1106 bp
RefSeq ORF:	810 bp
Locus ID:	66603
UniProt ID:	Q9CQQ4
Cytogenetics:	12 C1



[View online »](#)

Gene Summary:

This gene encodes one of the proteins found in the survival of motor neuron (SMN) complex, which consists of the SMN protein and several gemin proteins. The SMN complex is localized to a subnuclear compartment called gems (gemini of coiled bodies) and is required for assembly of spliceosomal small nuclear ribonucleoproteins (snRNP) and for pre-mRNA splicing. This protein interacts directly with the SMN protein and it is required for formation of the SMN complex. Disruption of this gene in mouse resulted in impaired snRNP assembly, and motor neuron degeneration. [provided by RefSeq, Sep 2015]