

## Product datasheet for **MR227475L3V**

### Tshr (NM\_001113404) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Tshr (NM_001113404) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tshr
Synonyms:	AI481368; hypothroid; hyt; pet
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113404
ORF Size:	735 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227475).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001113404.1</a> , <a href="#">NP_001106875.1</a>
RefSeq Size:	1346 bp
RefSeq ORF:	738 bp
Locus ID:	22095
Cytogenetics:	12 44.51 cM
Gene Summary:	Receptor for the thyroid-stimulating hormone (TSH) or thyrotropin. Also acts as a receptor for the heterodimeric glycoprotein hormone (GPHA2:GPHB5) or thyrostimulin. The activity of this receptor is mediated by G proteins which activate adenylate cyclase. Plays a central role in controlling thyroid cell metabolism.[UniProtKB/Swiss-Prot Function]



[View online »](#)