

Product datasheet for **MR227482A1V**

Mouse Fgf23 (NM_022657) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Mouse Fgf23 (NM_022657) AAV Particle
Tag:	Myc-DDK
Symbol:	Fgf23
Synonyms:	Fgf8b
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>MR227482 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTAGGGACCTGCCTTAGACTCCTGGTGGGCGCGCTCTGCACTGTCTGCAGCTTGGGCACTGCTAGAG
CCTATCCAGACACTTCCCCATTGCTTGGCTCCAAGTGGGAAGCCTGACCCACCTGTACACGGCTACAGC
CAGGACCAGCTATCACCTACAGATCCATAGGGATGGTCATGTAGATGGCACCCCCATCAGACCATCTAC
AGTGCCCTGATGATTACATCAGAGGACGCCGGCTCTGTGGTGATAACAGGAGCCATGACTCGAAGTTCC
TTTGTATGGATCTCCACGGCAACATTTTGGATCGCTTCACTTCAGCCAGAGAATTGCAAGTTCCGCCA
GTGGACGCTGGAGAATGGCTATGACGTCTACTTGTGCGAGAAGCATCACTACCTGGTGAGCCTGGGCCCG
GCCAAGCGCATCTTCCAGCCGGGCACCAACCCGCCCTTCTCCAGTTCTTGCGCGCAGGAACGAGG
TCCCGTGCTGCACTTCTACACTGTTGCGCCACGGCGCCACACGCGCAGCGCCGAGGACCCACCCGAGCG
CGACCCACTGAACGTGCTCAAGCCGCGGCCCGCGCCACGCCTGTGCCTGTATCCTGCTCTCGCGAGCTG
CCGAGCGCAGAGGAAGGTGGCCCCGAGCCAGCGATCCTCTGGGGGTGCTGCGCAGAGGCCGTGGAGATG
CTCGCGGGGGCGCGGGAGGCGCGGATAGGTGTCGCCCCTTTCCAGGTTTCGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR227482 protein sequence Red=Cloning site Green=Tags(s) MLGTCLRLLVGALCTVCSLGTARAYPDTSPLLGSNWGSLTHLYTATARTSYHLQIHRDGHVDGTPHQTIY SALMITSEDAGSVVITGAMTRRFLCMDLHGNI FGSLHFSPENCKFRQWTL ENG YD VYLSQKHHYLVSLGR AKRIFQPGTNPPFSQFLARRNEVPLLHFYTVRPRRHTRSAEDPPERDPLNVLKPRPRATPVVSCSREL PSAEEGGPAASDPLGVLRRGRGDARGGAGGADRCRPFPRFV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_022657
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_022657.3, NP_073148.1</u>
RefSeq Size:	2702 bp
RefSeq ORF:	756 bp
Locus ID:	64654
UniProt ID:	<u>Q9EPC2</u>
Cytogenetics:	6 F3
MW:	27.7 kDa