

Product datasheet for **MR227495L3V**

Atxn3 (NM_029705) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Atxn3 (NM_029705) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Atxn3
Synonyms:	2210008M02Rik; AI463012; AI647473; ataxin-3; ATX3; Mjd; MJD1; Sca3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029705
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227495).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029705.3 , NP_083981.2
RefSeq Size:	5376 bp
RefSeq ORF:	1068 bp
Locus ID:	110616
UniProt ID:	Q9CVD2
Cytogenetics:	12 E



[View online »](#)

Gene Summary:

Deubiquitinating enzyme involved in protein homeostasis maintenance, transcription, cytoskeleton regulation, myogenesis and degradation of misfolded chaperone substrates (By similarity). Binds long polyubiquitin chains and trims them, while it has weak or no activity against chains of 4 or less ubiquitins (By similarity). Involved in degradation of misfolded chaperone substrates via its interaction with STUB1/CHIP: recruited to monoubiquitinated STUB1/CHIP, and restricts the length of ubiquitin chain attached to STUB1/CHIP substrates and preventing further chain extension (PubMed:21855799). Interacts with key regulators of transcription and represses transcription: acts as a histone-binding protein that regulates transcription (By similarity). Regulates autophagy via the deubiquitination of 'Lys-402' of BECN1 leading to the stabilization of BECN1 (PubMed:28445460).[UniProtKB/Swiss-Prot Function]