

## Product datasheet for **MR227624L3V**

### Phospho1 (NM\_153104) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Phospho1 (NM_153104) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Phospho1
Synonyms:	D11Moh36; Phospo1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153104
ORF Size:	801 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR227624).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153104.3</a> , <a href="#">NP_694744.1</a>
RefSeq Size:	1910 bp
RefSeq ORF:	804 bp
Locus ID:	237928
UniProt ID:	<a href="#">Q8R2H9</a>
Cytogenetics:	11 59.01 cM
Gene Summary:	Phosphatase that has a high activity toward phosphoethanolamine (PEA) and phosphocholine (PCho). Involved in the generation of inorganic phosphate for bone mineralization.[UniProtKB/Swiss-Prot Function]



[View online »](#)