

Product datasheet for **MR227626**

Pfkfb3 (NM_001177753) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pfkfb3 (NM_001177753) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pfkfb3
Synonyms:	E330010H22Rik; iPFK-2; uPFK-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR227626 representing NM_001177753
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGCCGTTGGAAGTACCCAGAGCCGGTACAGAAGATCTGGGTGCCCGTCGATCACCGACCTCGTTGCCAGATCCTGTGGGCCGAAGCTGACCAACTCCCAACCGTGATTGTGATGGTGGGCTCCCAGCCGGGTAAGACTTACATCTCAAGAAGCTGACTCGCTACCTCAACTGGATAGGTGTTCCAACGAAAGTGTCAATGTGGGAGAGTATCGCGTGAGGCTGTGAAGCAGTACAGCTTTACAACCTCTCCGCCCTGACAATGAGGAAGCCATGAGAGTCCGAAACAGTGTGCCTTAGCTGCCCTGAGAGATGTCAAAGCTACCTGACAAAGGAAGGTGGACAGATTGCAGTTTTTGTGCCACCAATACAACCTAGAGAAAGAAGGCACATGATCCTTAATTTGCCAAAGAAAATGACTTCAAGGCATTTTTCATCGAGTCGGTCTGTGACGACCCTACTGTTGTGCCCTCCACATCATGGAAGTAAAACTCCAGCCCGATTACAAAGACTGCAACTCAGCAGAAGCCATGGACGACTTCATGAAGAGAATCAATTGCTATGAAGCCAGCTACCAGCCTTTGACCCTGATAAATGTGACAGGGACTTGTACTCATCAAAGTGATTGATGTTGGCCGGAGGTTCTTGGTGAACAGAGTGCAGGACCACATCCAGAGCCGAATTGTATACTACCTGATGAACATCCAGTGCAGCCCGTACCATCTACCTGTGTGCGCATGGCGAGAA TGAGTACAACCTTGAGGGCAAGATTGGGGGCGACTCAGGCTTGTCCAGCAGAGGCAAGAAGTTCGCCAATGCCCTGAGCAAGTTCGTAGAAGAACAGAATCTGAAGGACCTCCGTGTGTGGACCAGCCAGCTGAAAAGTACCATTAGACCCGGAGGCCCTGAGGCTCCCCTACGAGCAGTGAAGGCACTCAACGAGATCGATGCTGGTGTGTGAGGAACTAACCTATGAGGAGATCAGAGATACCTACCTGAGGAGTACGCACTGCGTGAACAGACAAGTACTACTATCGCTATCCACGGGAGATCCTACCAGGATCTTGTCCAACGCCTGGAGCCTGTGATCAGGAAATAGAGCGCAAGAGAATGTATTGGTCACTGTACCAGGCTGTTCTACGCTGCCTGATGCCTACTTCTCGACAAGAGTGCAGAGGAGATGCCATACCTGAAATGTCCGCTCCACACTGTCTGAAACTGACACCTGTGCATATGGGTGCCGGTGGAGTCCATTTACTTGAATGTAGAATCGGTGAGCACACCCGGGAGAGGTGAGAGGCTGTCAAAAACAACACTTTGCCAGTGTAGTCAGACCTTCTCTTACACCGAACTTGA CTTTCAGTCTGTGGAGTCCGCAAAACAGGATGCAAAAGAAGGGACCTAACCCGCTCATGAGACGCAATAGTGTACACCCACTAGCCAGTCCCGAGCCACCAAAAAGCCTCGCATCAACAGCTTTGAGGAGCGTGTGGCTCCACTTCTGCTGCCCTGCCAAGCTGCCTGCCCCGGAAGTGCCACGCAGCTGCCCGGACAACCTTTGCTAGGAAAGCCTGTCTACGAAGTGTCTGTACATTTTTTCAAAGTTTTCTCCTTAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR227626 representing NM_001177753
 Red=Cloning site Green=Tags(s)

MPLELTQSRVQKIWVPVDHRPSLPRSCGPKLTNSPTVIVMVGLPARGKTYISKKLTRYLNWIGVPTKVFNVGEYRREAVKQYSSYNFFRPDNEAMRVRKQCALAALRDVKSYLEKEGGQIAVFDATNTTRERRHMILNFAKENDFKAFFIESVDDPTVVASNIMEVKISSPDYKDCNSAEAMDDFMKRINCYEASYQLDPDKDRDL SLIKVIDVGRFLVNRVQDHIQSRIVVYLMNIHVQPRTIYLCRHGENEYNLQGKIGGDSGLSSRGKFFANALSKFVEEQNLKDLRVWTSQKSTIQTAELRLPYEQWKALNEIDAGVCEELTYEIRDTYPEEYALREQDKYYYYRYPTGESYQDLVQRLEPVIMELERQENLVICHQAVLRCLLAYFLDKSAEEMPYLKPLHTVLKLTTPVAYGCRVESIYLNVESVSTHRERSEAVKIQHFASVVRPSSYTELDFQSVESAKQDAKGNPLMRRNSVTPLASPEPTKKPRINSFEERVASTSAALPSCLPPEVPTQLPGQPLLGKACLRSVCHIFSKFSKY

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001177753

ORF Size: 1665 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177753.1](#), [NP_001171224.1](#)

RefSeq Size: 5206 bp


RefSeq ORF: 1668 bp

Locus ID: 170768

Cytogenetics: 2 A1

MW: 64.1 kDa

Product images:



Circular map for MR227626