

Product datasheet for **MR227632**

Pfkfb3 (NM_001177755) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pfkfb3 (NM_001177755) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pfkfb3
Synonyms:	E330010H22Rik; iPFK-2; uPFK-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR227632 representing NM_001177755
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCGTTGGAAGTACCCAGAGCCGGGTACAGAAGATCTGGGTGCCCGTCGATACCGACCTCGTTGC
 CCAGATCTGTGGGCCGAAGCTGACCAACTCCCAACCGTATTGTGATGGTGGGCTCCAGCCGGGG
 TAAGACTTACATCTCAAGAAGCTGACTCGCTACCTCACTGGATAGGTGTTCCAACGAAAGTGTCAAT
 GTGGGAGAGTATCGGCGTGAGGCTGTGAAGCAGTACAGCTTTACAACCTCTCCGCCCTGACAATGAGG
 AAGCCATGAGAGTCCGAAACAGTGTGCCTTAGCTGCCCTGAGAGATGTCAAAGCTACCTGACAAAGGA
 AGGTGGACAGATTGCAGTTTTTGTGCCACCAATACAACCTAGAGAAAGAAGGCACATGATCCTTAATTTT
 GCCAAAGAAAATGACTTCAAGGCATTTTTCATCGAGTCGGTCTGTGACGACCCTACTGTTGTGCCCTCCA
 ACATCATGGAAGTAAAACTCCAGCCCGATTACAAAGACTGCAACTCAGCAGAAGCCATGGACGACTT
 CATGAAGAGAATCAATTGCTATGAAGCCAGCTACCAGCCTTTGACCCTGATAAATGTGACAGGGACTTG
 TCACTCATCAAAGTGATTGATGTTGGCCGGAGGTTCTTGGTGAACAGAGTGCAGGACCACATCCAGAGCC
 GAATTGTATACTACCTGATGAACATCCAGTGCAGCCCGTACCATCTACCTGTGTCGGCATGGCGAGAA
 TGAGTACAACCTTGAGGGCAAGATTGGGGGCGACTCAGGCTTGTCCAGCAGAGGCAAGAAGTTCGCCAAT
 GCCCTGAGCAAGTTCGTAGAAGAACAGAATCTGAAGGACCTCCGTGTGTGGACCAGCCAGCTGAAAAGTA
 CCATTCAGACCGCGGAGGCCCTGAGGCTCCCCTACGAGCAGTGAAGGCACTCAACGAGATCGATGCTGG
 TGTGTGTGAGGAACTAACCTATGAGGAGATCAGAGTACCTACCTGAGGAGTACGCACTGCGTGAACAG
 GACAAGTACTACTATCGCTATCCACCGGAGAGTCCCTACCAGGATCTTGTCCAACGCCTGGAGCCTGTGA
 TCATGGAATTAGAGCGCCAAGAGAATGTATTGGTCACTGTACCAGGCTGTTCTACGCTGCCTACTAGC
 CTACTTCTCGACAAGAGTGCAGAGGAGATGCCATACCTGAAATGTCCGCTCCACACTGTCTGAAACTG
 ACACCTGTGGCATATGGGTGCCGGTGGAGTCCATTTACTTGAATGTAGAATCGGTGAGCACACCCGGG
 AGAGGTGAGAGGCTGTCAAAATACAACCTTTGCCAGTGTAGTCAGACCTTCTCTTACACCGAAGTGA
 CTTTCAGTCTGTGGAGTCCGCAAAACAGGATGCAAGAAGGGACCTAACCCGCTCATGAGACGCAATAGT
 GTCACCCCACTAGCCAGTCCCGAGCCACCAAAAAGCCTCGCATCAACAGCTTTGAGGAGCGTGTGGCTT
 CCATTCTGCTGCCCTGCCAAGCTGCCTGCCCCGGAAGTGCCACGCAGCTGCCCGGACAAAACATGAA
 GAGCTCCCGGAGCGGTGCCGACTTCCAGAAGCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR227632 representing NM_001177755
 Red=Cloning site Green=Tags(s)

MPLELTQSRVQKIWVPVDHRPSLPRSCGPKLTNSPTVIVMVGLPARGKTYISKKLTRYLNWIGVPTKVFN
 VGEYRREAVKQYSSYNFFRPDNEAMRVRKQCALAALRDVKSYLEKEGGQIAVFDATNTTRERRHMILNF
 AKENDFKAFFIESVDDPTVVASNIMEVKISSPDYKDCNSAEAMDDFMKRINCYEASYQLDPDKDRDL
 SLIKVIDVGRFLVNRVQDHIQSRIVYYLMIHVQPRTIYLCRHGENEYNLQGKIGGDSGLSSRGKFFAN
 ALSKFVEEQNLKDLRVWTSQKSTIQTAELRPLPYEQWKALNEIDAGVCEELTYEEIRDTYPEEYALREQ
 DKYYYYRYPTGESYQDLVQRLEPVIMELERQENLVICHQAVLRCLLAYFLDKSAEEMPYLKPLHTVLKL
 TPVAYGCRVESIYLNVESVSTHRERSEAVKIQHFASVVRPSSYTELDFQSVESAKQDAKGNPLMRNS
 VTPLASPEPTKPRINSFEERVASTSAALPSCLPPEVPTQLPGQNMKSSRSGADSSQKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/ja3487_a01.zip

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001177755

ORF Size: 1650 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177755.1](#), [NP_001171226.1](#)

RefSeq Size: 5130 bp
 RefSeq ORF: 1650 bp
 Locus ID: 170768
 Cytogenetics: 2 A1
 MW: 62.9 kDa

Product images:



Circular map for MR227632